

COOK®

MEDICAL

The Zenith® AAA Aortouni-Iliac Graft is continuing the COOK® tradition of offering superior products for improved patient care. Our commitment to education, research and development will keep you on the cutting edge of endovascular technology.

Zenith®

AAA AORTOUNI-ILIAC GRAFT



COOK®  
MEDICAL

[www.cookmedical.com](http://www.cookmedical.com)

COOK MEDICAL INCORPORATED  
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812 339-2235, Toll Free: 800 457-4500, Toll Free Fax: 800 554-8335

COOK (CANADA) INC.  
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905 640-7110, Toll Free: 800 668-0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.  
Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS  
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2008

AI-DM-ZAAA-AUI-0802-340-04EN

## TRANSCENDING AORTIC REPAIR

Success in treating acute and challenging cases requires fast and safe intervention. That is exactly what the Zenith® AAA Aortouni-iliac device is all about. Its modular one- or two- segment design allows for a perfect match even in difficult anatomies. The straightforward introducer system features a controlled release to help deploy the endovascular graft precisely where you want it – in a minimum of time.

### FIXATION

The Zenith AAA Aortouni-Iliac Graft offers **suprarenal fixation with anchoring barbs** which reduce the risk of migration and enhance graft-to-vessel sealing.

- Fixation barbs are positioned incrementally on the suprarenal stent to provide added stability.
- Specially designed z-stents also provide increased graft/vessel apposition to minimize endoleaks and graft migration.

### DESIGN

Distinct Zenith features provide you with the tools you need to **perform with confidence**.

- Independent, stainless steel z-stent design provides excellent flexibility and columnar strength.
- Zenith AAA Aortouni-Iliac Grafts are made with a lightweight yet strong shrink-resistant woven polyester with proven AAA treatment capabilities.

### DELIVERY

The Zenith AAA Aortouni-Iliac Graft with **H&L-B One-Shot™ Introduction System** offers a unique control method for endovascular graft deployment.

- Sequential deployment provides continuous control of the entire graft throughout the duration of the procedure.
- Enables precise positioning and the ability to readjust the endovascular graft before deploying the fixation barbs.

