




COMPONENT ORDERING INFORMATION

ZENITH FLEX AAA ENDOVASCULAR GRAFT with Z-Trak™ Introduction System

Bifurcated Main Body Graft




Global Product Number	Order Number	Main Body Diameter mm	Main Body Length mm	Introducer Sheath ID/OD
22 mm Diameter				
G48402	TFFB-22-82-ZT	22	82	18 Fr (6.0 mm)/7.1 mm
G48403	TFFB-22-96-ZT	22	96	18 Fr (6.0 mm)/7.1 mm
G48404	TFFB-22-111-ZT	22	111	18 Fr (6.0 mm)/7.1 mm
G48405	TFFB-22-125-ZT	22	125	18 Fr (6.0 mm)/7.1 mm
G48406	TFFB-22-140-ZT	22	140	18 Fr (6.0 mm)/7.1 mm
24 mm Diameter				
G48407	TFFB-24-82-ZT	24	82	18 Fr (6.0 mm)/7.1 mm
G48408	TFFB-24-96-ZT	24	96	18 Fr (6.0 mm)/7.1 mm
G48409	TFFB-24-111-ZT	24	111	18 Fr (6.0 mm)/7.1 mm
G48410	TFFB-24-125-ZT	24	125	18 Fr (6.0 mm)/7.1 mm
G48411	TFFB-24-140-ZT	24	140	18 Fr (6.0 mm)/7.1 mm
26 mm Diameter				
G48412	TFFB-26-82-ZT	26	82	18 Fr (6.0 mm)/7.1 mm
G48413	TFFB-26-96-ZT	26	96	18 Fr (6.0 mm)/7.1 mm
G48414	TFFB-26-111-ZT	26	111	18 Fr (6.0 mm)/7.1 mm
G48415	TFFB-26-125-ZT	26	125	18 Fr (6.0 mm)/7.1 mm
G48416	TFFB-26-140-ZT	26	140	18 Fr (6.0 mm)/7.1 mm
28 mm Diameter				
G48417	TFFB-28-82-ZT	28	82	20 Fr (6.7 mm)/7.7 mm
G48418	TFFB-28-96-ZT	28	96	20 Fr (6.7 mm)/7.7 mm
G48419	TFFB-28-111-ZT	28	111	20 Fr (6.7 mm)/7.7 mm
G48420	TFFB-28-125-ZT	28	125	20 Fr (6.7 mm)/7.7 mm
G48421	TFFB-28-140-ZT	28	140	20 Fr (6.7 mm)/7.7 mm
30 mm Diameter				
G48422	TFFB-30-82-ZT	30	82	20 Fr (6.7 mm)/7.7 mm
G48423	TFFB-30-96-ZT	30	96	20 Fr (6.7 mm)/7.7 mm
G48424	TFFB-30-111-ZT	30	111	20 Fr (6.7 mm)/7.7 mm
G48425	TFFB-30-125-ZT	30	125	20 Fr (6.7 mm)/7.7 mm
G48426	TFFB-30-140-ZT	30	140	20 Fr (6.7 mm)/7.7 mm
32 mm Diameter				
G48427	TFFB-32-82-ZT	32	82	20 Fr (6.7 mm)/7.7 mm
G48428	TFFB-32-96-ZT	32	96	20 Fr (6.7 mm)/7.7 mm
G48429	TFFB-32-111-ZT	32	111	20 Fr (6.7 mm)/7.7 mm
G48430	TFFB-32-125-ZT	32	125	20 Fr (6.7 mm)/7.7 mm
G48431	TFFB-32-140-ZT	32	140	20 Fr (6.7 mm)/7.7 mm
36 mm Diameter				
G48432	TFFB-36-95-ZT	36	95	22 Fr (7.3 mm)/8.5 mm
G48433	TFFB-36-113-ZT	36	113	22 Fr (7.3 mm)/8.5 mm
G48434	TFFB-36-131-ZT	36	131	22 Fr (7.3 mm)/8.5 mm
G48435	TFFB-36-149-ZT	36	149	22 Fr (7.3 mm)/8.5 mm

Sheath length for all main body graft diameters = 40 cm.

Sheath and valve length for all main body graft diameters = 51 cm.

ZENITH FLEX AAA ENDOVASCULAR GRAFT with Z-Trak Introduction System

Contralateral/Ipsilateral Iliac Legs







Sheath length for all iliac leg graft diameters = 71 cm.

Sheath and valve length for all iliac leg graft diameters = 82 cm.

Global Product Number	Order Number	Distal Iliac Leg Diameter mm	Working Length of Iliac Leg mm	Introducer Sheath ID/OD
8 mm Diameter				
G48436	TFLE-8-37-ZT	8	37	14 Fr (4.7 mm)/5.3 mm
G48437	TFLE-8-54-ZT	8	54	14 Fr (4.7 mm)/5.3 mm
G48438	TFLE-8-71-ZT	8	71	14 Fr (4.7 mm)/5.3 mm
G48439	TFLE-8-88-ZT	8	88	14 Fr (4.7 mm)/5.3 mm
G48440	TFLE-8-105-ZT	8	105	14 Fr (4.7 mm)/5.3 mm
G48441	TFLE-8-122-ZT	8	122	14 Fr (4.7 mm)/5.3 mm
10 mm Diameter				
G48442	TFLE-10-37-ZT	10	37	14 Fr (4.7 mm)/5.3 mm
G48443	TFLE-10-54-ZT	10	54	14 Fr (4.7 mm)/5.3 mm
G48444	TFLE-10-71-ZT	10	71	14 Fr (4.7 mm)/5.3 mm
G48445	TFLE-10-88-ZT	10	88	14 Fr (4.7 mm)/5.3 mm
G48446	TFLE-10-105-ZT	10	105	14 Fr (4.7 mm)/5.3 mm
G48447	TFLE-10-122-ZT	10	122	14 Fr (4.7 mm)/5.3 mm
12 mm Diameter				
G48448	TFLE-12-39-ZT	12	39	16 Fr (5.3 mm)/6.0 mm
G48449	TFLE-12-56-ZT	12	56	16 Fr (5.3 mm)/6.0 mm
G48450	TFLE-12-73-ZT	12	73	16 Fr (5.3 mm)/6.0 mm
G48451	TFLE-12-90-ZT	12	90	16 Fr (5.3 mm)/6.0 mm
G48452	TFLE-12-107-ZT	12	107	16 Fr (5.3 mm)/6.0 mm
G48453	TFLE-12-124-ZT	12	124	16 Fr (5.3 mm)/6.0 mm
14 mm Diameter				
G48454	TFLE-14-39-ZT	14	39	16 Fr (5.3 mm)/6.0 mm
G48455	TFLE-14-56-ZT	14	56	16 Fr (5.3 mm)/6.0 mm
G48456	TFLE-14-73-ZT	14	73	16 Fr (5.3 mm)/6.0 mm
G48457	TFLE-14-90-ZT	14	90	16 Fr (5.3 mm)/6.0 mm
G48458	TFLE-14-107-ZT	14	107	16 Fr (5.3 mm)/6.0 mm
G48459	TFLE-14-124-ZT	14	124	16 Fr (5.3 mm)/6.0 mm
16 mm Diameter				
G48460	TFLE-16-39-ZT	16	39	16 Fr (5.3 mm)/6.0 mm
G48461	TFLE-16-56-ZT	16	56	16 Fr (5.3 mm)/6.0 mm
G48462	TFLE-16-73-ZT	16	73	16 Fr (5.3 mm)/6.0 mm
G48463	TFLE-16-90-ZT	16	90	16 Fr (5.3 mm)/6.0 mm
18 mm Diameter				
G48464	TFLE-18-39-ZT	18	39	16 Fr (5.3 mm)/6.0 mm
G48465	TFLE-18-56-ZT	18	56	16 Fr (5.3 mm)/6.0 mm
G48466	TFLE-18-73-ZT	18	73	16 Fr (5.3 mm)/6.0 mm
G48467	TFLE-18-90-ZT	18	90	16 Fr (5.3 mm)/6.0 mm
20 mm Diameter				
G48468	TFLE-20-39-ZT	20	39	16 Fr (5.3 mm)/6.0 mm
G48469	TFLE-20-56-ZT	20	56	16 Fr (5.3 mm)/6.0 mm
G48470	TFLE-20-73-ZT	20	73	16 Fr (5.3 mm)/6.0 mm
G48471	TFLE-20-90-ZT	20	90	16 Fr (5.3 mm)/6.0 mm
22 mm Diameter				
G48472	TFLE-22-39-ZT	22	39	16 Fr (5.3 mm)/6.0 mm
G48473	TFLE-22-56-ZT	22	56	16 Fr (5.3 mm)/6.0 mm
G48474	TFLE-22-73-ZT	22	73	16 Fr (5.3 mm)/6.0 mm
G48475	TFLE-22-90-ZT	22	90	16 Fr (5.3 mm)/6.0 mm
24 mm Diameter				
G48476	TFLE-24-39-ZT	24	39	16 Fr (5.3 mm)/6.0 mm
G48477	TFLE-24-56-ZT	24	56	16 Fr (5.3 mm)/6.0 mm
G48478	TFLE-24-73-ZT	24	73	16 Fr (5.3 mm)/6.0 mm
G48479	TFLE-24-90-ZT	24	90	16 Fr (5.3 mm)/6.0 mm

ZENITH FLEX AAA ENDOVASCULAR GRAFT

Ancillary Components



Global Product Number	Order Number	Proximal Diameter mm	Distal Diameter mm	Length mm	Introducer Sheath ID/OD
Converters					
G12496	ESC-24-12-80	24	12	80	18 Fr (6.0 mm)/6.9 mm
G12483	ESC-28-12-80	28	12	80	20 Fr (6.7 mm)/7.7 mm
G12484	ESC-32-12-80	32	12	80	20 Fr (6.7 mm)/7.7 mm
G43355	ESC-36-12-82	36	12	82	20 Fr (6.7 mm)/7.7 mm

Global Product Number	Order Number	Diameter mm	Length mm	Introducer Sheath ID/OD
Main Body Extensions				
G32179	ESBE-22-39	22	39	18 Fr (6.0 mm)/6.9 mm
G32180	ESBE-24-39	24	39	18 Fr (6.0 mm)/6.9 mm
G32181	ESBE-26-39	26	39	18 Fr (6.0 mm)/6.9 mm
G32182	ESBE-28-39	28	39	20 Fr (6.7 mm)/7.7 mm
G32183	ESBE-30-39	30	39	20 Fr (6.7 mm)/7.7 mm
G32184	ESBE-32-39	32	39	20 Fr (6.7 mm)/7.7 mm
G43360	ESBE-36-50	36	50	20 Fr (6.7 mm)/7.7 mm

Global Product Number	Order Number	Diameter mm	Length mm	Introducer Sheath ID/OD
Main Body Extensions				
G42783	ESBE-22-58	22	58	18 Fr (6.0 mm)/6.9 mm
G42784	ESBE-24-58	24	58	18 Fr (6.0 mm)/6.9 mm
G42785	ESBE-26-58	26	58	18 Fr (6.0 mm)/6.9 mm
G42786	ESBE-28-58	28	58	20 Fr (6.7 mm)/7.7 mm
G42787	ESBE-30-58	30	58	20 Fr (6.7 mm)/7.7 mm
G42788	ESBE-32-58	32	58	20 Fr (6.7 mm)/7.7 mm
G43361	ESBE-36-73	36	73	20 Fr (6.7 mm)/7.7 mm

Global Product Number	Order Number	Diameter mm	Length mm	Introducer Sheath ID/OD
Iliac Plugs				
G39468	ZIP-14	14	30	14 Fr (4.7 mm)/5.3 mm
G21365	ZIP-16	16	30	14 Fr (4.7 mm)/5.3 mm
G39470	ZIP-20	20	30	16 Fr (5.3 mm)/6.0 mm
G21366	ZIP-24	24	30	16 Fr (5.3 mm)/6.0 mm



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812 339-2235, Toll Free: 800 457-4500, Toll Free Fax: 800 554-8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905 640-7110, Toll Free: 800 668-0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2009

AI-DM-FXZTCOI-0905-508-02EN

AORTIC
INTERVENTION

CARDIOLOGY

CRITICAL
CARE

ENDOSCOPY

INTERVENTIONAL
RADIOLOGY

PERIPHERAL
INTERVENTION

SURGERY

UROLOGY

WOMEN'S
HEALTH