




Zenith Flex[®]

AAA ENDOVASCULAR GRAFT

PRODUCT BULLETIN

Flex Legs—Design and Delivery

The redesigned Zenith Flex Legs are now more kink-resistant and come in a wider variety of sizes. The look and feel have changed as well with the ergonomic Z-Trak™ Introduction System, featuring a new Captor[®] valve and a longer sheath.



DESIGN



Flex Design

Shorter stents and increased gap length give the Flex legs 24%-35% decrease in kink radius.



Options

More sizes have also been added to expand your treatment options.

DELIVERY



Z-Trak Introduction System

Redesigned Captor valve is adjustable and inhibits blood reflux. Increased sheath length complements the main body introduction system, easing delivery of the Flex leg graft.



www.cookmedical.com

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.
Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2008 AI-BOUSM-FXLGPB-EN-200806