




# Zenith® Renu™

AAA ANCILLARY GRAFT

## COMPONENT ORDERING INFORMATION

### ZENITH RENU AAA ANCILLARY GRAFT

#### Converter




Global Product Number	Order Number	Proximal Diameter mm	Distal Diameter mm	Introducer Sheath French Size	Proximal Length mm	Distal Length mm	Total Length mm
<b>22 mm Diameter</b>							
G36147	AX1-1-22-108	22	12	18	37	71	108
G36148	AX1-1-22-125	22	12	18	37	88	125
<b>24 mm Diameter</b>							
G36153	AX1-1-24-108	24	12	18	37	71	108
G36154	AX1-1-24-125	24	12	18	37	88	125
<b>26 mm Diameter</b>							
G36159	AX1-1-26-108	26	12	18	37	71	108
G36160	AX1-1-26-125	26	12	18	37	88	125
<b>28 mm Diameter</b>							
G36165	AX1-1-28-108	28	12	20	37	71	108
G36166	AX1-1-28-125	28	12	20	37	88	125
<b>30 mm Diameter</b>							
G36171	AX1-1-30-108	30	12	20	37	71	108
G36172	AX1-1-30-125	30	12	20	37	88	125
<b>32 mm Diameter</b>							
G36177	AX1-1-32-108	32	12	20	37	71	108
G36178	AX1-1-32-125	32	12	20	37	88	125
<b>36 mm Diameter</b>							
G43370	AX1-1-36-116	36	12	22	52	64	116

**Introducer Sheaths** — (22 - 26 mm components) = 18 French = 7.1 mm OD  
 (28 - 32 mm components) = 20 French = 7.7 mm OD  
 (36 mm components) = 22 French = 8.5 mm OD

# ZENITH RENU AAA ANCILLARY GRAFT

## Main Body Extension



Global Product Number	Order Number	Diameter mm	Length mm	Introducer Sheath French Size
<b>22 mm Diameter</b>				
G36128	RX1-22-43	22	43	18
G36129	RX1-22-62	22	62	18
<b>24 mm Diameter</b>				
G36131	RX1-24-43	24	43	18
G36132	RX1-24-62	24	62	18
<b>26 mm Diameter</b>				
G36134	RX1-26-43	26	43	18
G36135	RX1-26-62	26	62	18
<b>28 mm Diameter</b>				
G36137	RX1-28-43	28	43	20
G36138	RX1-28-62	28	62	20
<b>30 mm Diameter</b>				
G36140	RX1-30-43	30	43	20
G36141	RX1-30-62	30	62	20
<b>32 mm Diameter</b>				
G36143	RX1-32-43	32	43	20
G36144	RX1-32-62	32	62	20
<b>36 mm Diameter</b>				
G43362	RX1-36-54	36	54	22

Introducer Sheaths { (22 - 26 mm components) = 18 French = 7.1 mm OD  
 (28 - 32 mm components) = 20 French = 7.7 mm OD  
 (36 mm components) = 22 French = 8.5 mm OD



[www.cookmedical.com](http://www.cookmedical.com)

### COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.

Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

### COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA

Phone: 905.640.7110, Toll Free: 800.668.0300

### WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains

Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

### WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2008

AI-DM-RCOI-0802-451-01EN